



ELEVATE YOUR CLINIC TO THE HIGHEST STANDARD OF
EXCELLENCE WITH HTM EQUIPMENT.
NEW TECHNOLOGIES, EFFICIENCY, DESIGN AND
PROFESSIONAL PERFORMANCE.

(+55) 19 3808-7741
www.htmeletronica.com.br
comex@htmeletronica.com.br
comercial@htmeletronica.com.br

Av. Rio Nilo, 209 - Jardim Figueira
CEP: 13904-380 - Amparo - São Paulo, Brasil



/HTMELETR



/htmeletronica



/htmeletronica

We reserve the right to correct any graphical errors. Merely
illustrative images. Internal distribution.





HTM

You better
every day

We are a reference industry in the development and manufacture of electronic equipment aimed at the aesthetics and physiotherapy segment in Brazil. With more than 20 years in the market, our differentials are based on high technology, exclusive design and great efficiency of our products.

All equipment is designed and developed in its own factory, with 100% national technology!

Our modern industrial structure has 5 thousand m², built under the strictest standards of production and environmental responsibility. We are constantly expanding, always recognizing and valuing the ethical commitment and well-being of our employees.



Innovating in technology and creating high quality equipment is our mission. In this path of overcoming, we continue to support and value our dealers, customers and partners with the commitment to offer the market the best and most effective products, integrating scientific knowledge with practicality.

Today, HTM Eletrônica is a leader in innovation, results and safety, parameters that will undoubtedly follow with us in the future.



Compact radiofrequency for body, facial and intimate treatments;

Quick reach of therapeutic temperature / Guarantee of maintenance of therapeutic temperature even in large areas;

Multifrequency Technology : 0.6MHz, 1.2MHz and 2.4MHz;

Monopolar, bipolar and multipolar application / Resistive and capacitive tips;

Innovative Pro Performance System, thought about the user experience.

SKIN FLACCIDITY;
FACIAL AGING AND WRINKLES;
ACNE SCARS;
STRETCH MARKS;
CELLULITIS;
LOCALIZED ADIPOSITY;
VULGAR ACNE;
INTIMATE AESTHETICS.

ANVISA - 80212480015

Two simultaneous and independent outputs;

Configurable cooling temperature down to -11°C (± 20%);

Exclusive iProtect system: guaranteed security and performance;

INFO menu with usage and care tutorial for better use of all equipment resources;

Intelligent programming with targeted protocols according to the treatment region and male and female option;

touch display screen in the cabinet and applicators;
Applicator cleaning system;

Exclusive i-check equipment integrity system;

Sophisticated design/Full range.

LOCALIZED ADIPOSITY

ANVISA - 80212480025
*All applicators are optional.



EFFECT

Radio frequency

Compact
high-performance
radiofrequency.

Conventional and contrast
cryolipolysis

Reduces up to 42% of
LOCALIZED ADIPOSITY in
just 1 session.

BEAUTY SHAPE DUO

LÍMINE

Radio frequency

Excellence and high
performance in
Radiofrequency
treatments.

Carboxytherapy

Intelligent system of
percutaneous and
transcutaneous
carboxytherapy.

PLURIA



Sophisticated design / Exclusive applicators for body, facial and intimate treatments;

Multifrequency Technology : 0.6MHz, 1.2MHz and 2.4MHz;

Multi Layer Treatment : automatic frequency variation, reaching different tissue depths;

Monopolar, bipolar and multipolar application / Resistive and capacitive tips;

Different emission modes: continuous and pulsed;
Quick reach of therapeutic temperature / Guarantee of maintenance of therapeutic temperature even in large areas;

Innovative Pro Performance System: thought about the user experience.

Sagging of the SKIN; FACIAL AGING AND WRINKLES;
ACNE SCARS;
STRETCH MARKS;
CELLULITIS;
LOCALIZED ADIPOSITY;
VULGAR ACNE; INTIMATE AESTHETICS.

ANVISA - 80212480015

Progressive gas heating and Comfort System;

Exclusive Carbo System Cupping for non-invasive carboxytherapy;

Pre-programmed protocols and possibility to save new ones;

Flow up to 400ml/min.

CELLULITIS;
ADIPOSITY;
SKIN FLACCIDITY.

ANVISA - 80212480020



LIGHT PULSE

Intense pulsed
light

Versatility and
high power in
light treatments.

Lamp with a guarantee of 1 million* activations;

Spot mode and scan mode*;

Electronic cooling "Crycooling";

New screen layout and interactive navigation, with the possibility of saving protocols;

Suggestion of treatment areas;

Epilation, rejuvenation, acne, wrinkles, blemishes, stretch marks and vascular therapy.

ACNE;
ACTINOSENESCENCE;
CHRONOSENESCENCE;
EPILATING;
STRETCH MARKS;
ROSACEA AND TELANGIECTASIA.

ANVISA 80212480021

*Applicator lamp is warranted for 1 million shots or 18 months, whichever comes first.

*In epilation and rejuvenation treatments.

VIBRIA MAXX



Exclusive



MODO
DUO

Ultrasound and
therapeutic
currents

High-performance
platform for
simultaneous
therapies and
low-frequency.

Possibility of performing individual therapies, combined therapies and simultaneous therapies - Exclusive DUO MODE;

1MHz and 3MHz body ultrasound with higher output power: 63W and 21cm² ERA - FloatingHead ® Technology;

Low Frequency Ultrasound with 40kHz, ERA of 10cm² and 30W of power;

Facial ultrasound - First 5MHz facial applicator in Brazil with 9W power and 3cm² ERA;

8 current channels: HighForce (1000Hz), Russian (2500Hz), MediumForce (4000Hz), SoftForce (8000Hz), Drainage and Electrolipolysis;

Exclusive HTM Smart ID System;

Pro Performance System;

INFO menu with usage and care tutorial for better use of all equipment resources;

4 channels for currents; Iontophoresis and MENS;

More than 100 pre-programmed protocols;

Complete line of accessories;

Sophisticated design with touch display 8" screen.

LOCALIZED ADIPOSITY;
CELLULITIS AND WATER RETENTION;
REJUVENATION AND STRETCH MARKS;
CUTANEOUS PERMEATION OF ASSETS;
REHABILITATION AND ANALGESIA;
TISSUE REPAIR;
AID IN MUSCLE STRENGTHENING (MUSCULAR Sagging).

ANVISA - 80212480023

Ultrasound and therapeutic currents;

Compact hybrid platform for simultaneous therapies;

1MHz and 3MHz body ultrasound with higher output power: 63W and ERA of 21cm²;

FloatingHead® Technology;

Possibility of performing individual therapies, combined therapies and simultaneous therapies
-Exclusive DUO MODE.

LOCALIZED ADIPOSITY;
CELLULITIS AND WATER RETENTION;
REJUVENATION AND STRETCH MARKS;
CUTANEOUS PERMEATION OF ASSETS;
REHABILITATION AND ANALGESIA;
TISSUE REPAIR;
MUSCLE STRENGTHENING.

ANVISA - 80212480023



HÍBRIDI

Ultrasound and
therapeutic currents

Compact hybrid platform
for simultaneous therapies.



MODO
DUO

Body ultrasound: 1 and 3 MHz, 24W of power and ERA of 8cm²;

Facial ultrasound: 5MHz, 9W of power and ERA of 3cm² (optional);


Currents: TENS, FES, HighForce (1000Hz), Microgalvanic, MENS, Interferential 2kHz and 4kHz, Calvanic, Russian Current (2500Hz), MediumForce (4000Hz), SoftForce (8000Hz) and Lipolysis;

2 channels for currents;

150 pre-programmed protocols.

ANALGESIA AND REHABILITATION;
CELLULITIS AND LOCALIZED ADIPOSITY;
STRETCH MARKS AND TISSUE REPAIR;
STRENGTHENING MUSCULAR;
REJUVENATION;
WATER RETENTION;
ASSET PERMEATION.

ANVISA - 8021280001



SONIC COMPACT MAXX

new ULTRAFOCUS

Body Focused Ultrasound,
Body Piezoelectric Shock
Waves, and Stationary
and Dynamic
Microfocused Facial
Ultrasound

Platform for Face Lifting
and Body Lipomodeling.

New model of body applicator with 4 tips options:
10mm, 15mm, 20mm and 30mm in depth;

New model of facial applicator with 3 tips options:
1.5mm, 3.0mm and 4.5mm in depth;

Sensor for decoupling the applicators on the skin
with an audible alert, which minimizes operational
errors in the application;

Automatic parameter adjustment according to the
selection of the tip depth and application
technique;

Stationary and dynamic application technique for
Facial Microfocalized Ultrasound;

Body piezoelectric shock waves with automatic
parameter adjustment based on the definition of
the treatment area;

Interactive navigation and intelligent selection of
treatment regions;

INFO menu and QR Code with access to use and
care tutorial for better use of all equipment
resources;

Pro Performance System;

8 inch Touch Screen Display;

Full Range;

Exclusive design.

LOCALIZED ADIPOSITY;
CELLULITIS AND STRETCH MARKS;
FACIAL AGING;
SKIN FLACCIDITY.

ANVISA - 80212480022





TECARE

Tecartherapy and radiofrequency

Challenge your senses with the new generation of Tecartherapy.



New Tecartherapy and Radiofrequency system;

Revolutionary "hands free" technique": transforming technology into a new sensory experience;

Exclusive Tecartherapy applicators specially developed for the perfect fit of the hands: ideal for the "hands on" technique";

Treatments in the aesthetic and rehabilitation area;

touch display 5" screen;

INFO menu with usage and care tutorial for better use of all equipment resources;

Intelligent programming with targeted protocols according to the therapeutic modality and treatment region;

Innovative Pro Performance System, thought about the user experience;

Exclusive applicator for intimate technique with remote activation and intensity adjustment;

Multifrequency Technology: 0.6MHz, 1.2MHz and 2.4MHz;

Multi Layer Treatment: automatic frequency variation, reaching different tissue depths;

Different emission modes: continuous and pulsed;

Different options of capacitive and resistive tips for Radiofrequency: monopolar, bipolar, tetrapolar, hexapolar and decapolar;

Tecartherapy tips of different sizes: 30,50,70 mm;

New plug for easy fitting and high durability;

Exclusive design;

Full range.

SAGGING OF THE SKIN;
FACIAL AGING AND WRINKLES;
ACNE SCARS;
STRETCH MARKS;
CELLULITIS;
LOCALIZED ADIPOSITY;

ACNE VULGAR;
INTIMATE
AESTHETICS;
INCONTINENCE
URINARY;
ANALGESIA;
REHABILITATION.

ANVISA - 80212480015

Allows you to choose between treatments: Vacuum LED, Vacuum Therapy, Phototherapy, Pressotherapy, Cupping Therapy, Diamond Peeling and Buttock Pump*;

Vacuum LED applicator containing: 3 green LEDs (300mW), 3 red LEDs (300mW) and 3 infrared LEDs (300mW), in continuous and pulsed modes;

It presents the possibility of individual or combined use of the LEDs, with or without the use of vacuum;

New key for Cupping Therapy, with a special and selective suction system for applying suction cups;

New digital pressure sensor with variable negative pressure adjustment from -1.5kPa to -80kPa**;

Pressotherapy with 2 output channels, with different modes of operation, insufflation and insufflation direction, being CHANNEL A for legs or arms and CHANNEL B for the abdomen;

Offers the best body, facial and capillary vacuum therapy applicators*;

INFO menu and QR Code with access to use and care tutorial for better use of all equipment resources;

More than 90 pre-programmed protocols and possibility to record new ones;

Language selection (Portuguese, English and Spanish);

Exclusive design with new digital pressure sensor;

5" screen touch display and medical standard plug for easy fitting and high durability;

Full range.

LOCALIZED ADIPOSITY AND CELLULITIS;
LYMPHATIC DRAINAGE;
FIBROSIS;
ANALGESIA;
RELAXING MASSAGE;
VENOUS STASIS;
LYMPHEDEMA AND EDEMA;
ACNE, STRETCH MARKS AND FACIAL AGING;
ALOPECIA AND MELASMA;
TISSUE REPAIR;
PHOTODYNAMIC THERAPY.

Anvisa - 80212480024

*Some accessories for performing these techniques are optional;

**Equivalent to -600mmHg.

Vacuum therapy,
cupping therapy,
pressotherapy and
phototherapy

Advanced
multifunctional
vacuum therapy ,
cupping therapy,
pressotherapy and
phototherapy
technology.

BEAUTY DERMO MAXX



SONIC COMPACT 1-3 MHZ

1 and 3 MHz, 21W of power and ERA of 7cm²;

24 pre-programmed protocols;

Intelligent programming of application time calculation (3MHz);

Connection for combination therapy.

CELLULITIS;
LOCALIZED ADIPOSITY;
WATER RETENTION;
TISSUE REPAIR;
ANALGESIA;
PHONOPHORESIS;
REHABILITATION.

ANVISA - 80212480001



Microcontrolled digital
ultrasound

High-tech portable
ultrasound.

Currents

Complete multiplatform
electrostimulation.



Microcontrolled digital
ultrasound

High technology in
portable ultrasound.

SONIC COMPACT 3MHZ

3 MHz, 24W of power and ERA of 8cm²;

Intelligent programming of application time calculation;

Automatic configuration of the ultrasonic intensity as a function of the adipose layer;

Connection for combined therapy;

11 pre-programmed protocols.

CELLULITIS;
LOCALIZED ADIPOSITY;
PHONOPHORESIS;
WATER RETENTION;
TISSUE REPAIR.

ANVISA - 80212480001

STIMULUS-R

TENS, FES, HighForce (1000Hz), Russian (2500Hz),
MediumForce (4000Hz), SoftForce (8000Hz),
Drainage, Lipolysis, Galvanic, Microgalvanic and
MENS;

10 channels for biphasic currents and 2 channels for
polarized currents;

46 pre-programmed protocols;

Easy access front exits;

Durable medical standard plug;

Current for electroacupuncture.

CELLULITIS;
LOCALIZED ADIPOSITY;
MUSCLE FLAVORITY;
CUTANEOUS PERMEATION OF ASSETS;
WATER RETENTION;
TISSUE REPAIR;
PHYSICAL REHABILITATION;
ANALGESIA.

ANVISA - 80212480011



Allows you to choose between treatments: Vacuum Therapy, Cupping Therapy, Diamond Peeling and Buttock Pump *;

New key for Cupping Therapy;

New digital pressure sensor;

Vacuum operating mode: Continuous and Pulsed;

It offers the best vacuum therapy applicators : body, facial and hair*;

More than 20 pre-programmed protocols;

QR Code with access to use and care tutorial for better use of all equipment resources.

LOCALIZED ADIPOSITY AND CELLULITIS;
LYMPHATIC DRAINAGE;
DECREASE OF FIBROSIS;
FACIAL AGING;
ANALGESIA AND RELAXING MASSAGE;
ACNE AND STRETCH MARKS.

ANVISA - 80212480024

*Some accessories for performing these techniques are optional.

TENS, FES, HighForce (1000Hz), Russian (2500Hz),
Lipolysis, Galvanic, Microgalvanic and MENS;

2 channels for currents;

28 pre-programmed body and facial protocols;

Top exits for easy access;

High resolution and contrast color OLED display.

MUSCULAR STRENGTHENING;
ANALGESIA AND FUNCTIONAL REHABILITATION;
DRAINAGE OF EDEMA;
CUTANEOUS PERMEATION OF ASSETS;
LOCALIZED ADIPOSITY AND CELLULITIS;
TISSUE REPAIR.

ANVISA - 80212480027



new BEAUTY DERMO

Vacuum therapy and
cupping therapy

Multifunctional vacuum
therapy and cupping
therapy platform for facial,
body and hair use.

Currents

Innovative portable
electrostimulator for
aesthetics.

STIM CARE



Look for a clinic with
HTM equipment. Quality
and safety in a brand I
trust!



FLUENCE MAXX is a Phototherapy equipment with the highest LED and LASER technology. It presents a wide range of resources for multitherapy and exclusive accessories, with 27 options: Laser Pens with higher power, LED Cluster (circular and linear), LED Cluster + Laser, Cluster LED MAXX, Cluster LED + Laser MAXX, Face LED Phototherapy Mask, Face LED Phototherapy Mask + Neck and Bracelet for ILIB Technique.

FLUENCE MAXX

LED and Laser

High technology in
LED and LASER
multitherapies.

Exclusive design;

5" touch screen display;

28 accessory options;

INFO menu with tutorial for to use and care for better use of all the equipment's resources;

Permit the individual choice or union of different light sources, coherent and non-coherent, increasing the spectrum of light spectrum and the range of treatments;

More than 1,700 * pre-programmed protocols and possibility to record new (* specific for each applicator, depending on the quantity of accessories purchased);

Toposcopic tip in the Laser Pens allowing to locate points for acupuncture;

Exclusive HTM Smart ID system for 904nm pens;

Medical standard plug of easy fitting and high durability;

Pen and cluster shutter button;

Exclusive bracelet for transcutaneous ILIB technique (optional);

LASER STOP and emergency button;

Laser / LED TEST device;

Full range.

ACNE;
ACUPUNCTURE;
ALOPECIA;
ANALGESIA;
CELLULITIS;
FAT LOCATED;
HYDRATION;
LYMPHATIC DRAINAGE;

LYMPHEDEMA;
MELASMA;
PHOTODYNAMIC
THERAPY;
PHYSICAL
REHABILITATION;
REJUVENATION;
TISSUE AND
INTRAVAGINAL REPAIR;
TRANSCUTANEOUS ILIB;
REJUVENATION.

ANVISA - 80212480028



Individual or associated therapies;

Option to include Laser pen, LED cluster and LED+Laser cluster;

More than 1,600 protocols*;

Exclusive HTM Smart ID system**.

ACNE;
ACUPUNCTURE - TRANSCUTANEOUS ILIB;
ALOPECIA;
ANALGESIA AND PHYSICAL REHABILITATION;
LOCALIZED ADIPOSITY AND CELLULITIS;
HYDRATION;
LYMPHEDEMA AND LYMPHATIC DRAINAGE;
MELASMA AND REJUVENATION;
TISSUE REPAIR;
PHOTODYNAMIC THERAPY.

ANVISA - 80212480028

*Depending on the number of accessories purchased.
**For 904nm pens.

Cluster with exclusive linear distribution;

Exclusive trichology suction cups for hair and brow applications;

Durable medical standard plug;

15 pre-programmed protocols.

ALOPECIA;
HAIR TREATMENTS;
SCALP PATHOLOGIES.

ANVISA - 80212480026

FLUENCE



LED and Laser

Innovative LED and Laser technology.



Phototherapy, vacuum therapy and high frequency

Platform for hair treatment.

STIM HAIR

ÁCRUS

Epilation laser with Triple Wave applicator (Alexandrite - 755nm, Diode-808nm, ND- Yag -1064nm)

Advanced Triple Wave LASER technology for epilation of all phototypes.



Meets all skin phototypes, including tanned ones;

Fast payback on investment;

Adaptation of wavelengths for all hair thicknesses and densities;

Smart Spot Cooling - MAXX Cryocooling technology down to -10 degrees - which guarantees greater comfort and safety for the patient;

Flexible treatment mode: Sweeping, Sequential and Punctual, ensuring application speed;

Pre-programmed parameters according to phototype, hair density and thickness or user adjustment option;

Equipment enabled with all the features of the exclusive HTM Match application;

Touch display 10.4" screen with greater precision to the touch and high durability;

Pro-Performance System with interactive navigation and selection of treatment regions;

INFO menu and QR Code with access to use and care tutorial for better use of all equipment resources;

Exclusive i-check equipment integrity system;

Activation mode: pedal (automatic or manual) or applicator (manual);

18-month warranty including unlimited shots* optical set;

Low-cost maintenance: no consumables;

Compact and exclusive design;

It comes with a carrying case that adapts to support the equipment and adapts to the environment.

EPILATION.

*Provided that the instructions described in the instruction manual regarding the technique of application, use, maintenance and cleaning are followed.

ANVISA - 80212480029

5 pre-programmed protocols;
Durable medical standard plug.

SKIN CLEANSING;
CELLULAR RENEWAL;
FACE LIFTING;
PEELING;
FACIAL DRAINAGE.

ANVISA - 8021280026

STIMULUS FACE CLEAN



Vacuum therapy ,
ultrasonic peeling and
high frequency

Skin cleansing platform.

STIMULUS FACE MAXX

HighForce (1000Hz), Russian (2500Hz), Galvanic,
Microgalvanic and MENS;

4 channels for biphasic currents and 2 channels for
polarized currents;

87 pre-programmed protocols;

Durable medical standard plug.

STRENGTHENING AND MUSCLE TONING;
FACIAL DRAINAGE AND PEELING;
SKIN CLEANING AND ACNE TREATMENT;
FACIAL LIFTING AND SKIN PERMEATION OF
ASSETS;
HYPERCHROMIAS AND CELL RENEWAL;
PHOTODYNAMIC THERAPY.

ANVISA - 8021280026

Innovative exclusive platform for skin cleaning;

HighForce (1000Hz), Russian (2500Hz), Galvanic,
Microgalvanic and MENS;


2 channels for biphasic and polarized currents;

85 pre-programmed protocols;

Durable medical standard plug.

MUSCLE STRENGTHENING AND TONING;
FACE DRAINAGE;
SKIN CLEANING AND ACNE TREATMENT;
FACIAL LIFTING AND SKIN PERMEATION OF
ASSETS;
PEELING AND TISSUE REPAIR;
PHOTODYNAMIC THERAPY.


ANVISA - 8021280026



Currents, facial vacuum
therapy, high frequency
and phototherapy

Face multiplatform.

STIMULUS FACE



New cabinet and applicator pen with unique design;

Automatic bivolt (1st and only in Brazil);

Long-range connection cable: offers greater freedom in handling;

Electrode holder (optional);

Universal fit for use with different electrodes;

Includes carrying bag.

SKIN TREATMENTS;
HEALING AND TISSUE REPAIR;
TREATMENTS IN PODOLOGY.

ANVISA - 80212480003



BEAUTY FACE

High frequency

The evolution in portable
High Frequency generator.

Revolutionary design, which combines beauty and practicality in operation;

It operates on two mains voltages - bivolt
(110V-/220V-);

Digital control of the water level of the reservoir,
against overflow;

Intelligent heating system that monitors the
temperature of the boiler providing more safety
and agility at the beginning of the operation;

Stem with 360° rotatable emitting tip and silicone
grip area for greater operator comfort;

Microcontroller with digital parameter controls;

Compartment for inserting flavorings.

ACNE;
ALOPECIA;
AROMATHERAPY;
HYDRATION.

ANVISA - 80212489004



BEAUTY STEAM new

Water steam, ozone
steam*, aromatherapy,
photobiomodulation and
chromotherapy

Skin hygiene,
aromatherapy and light
therapies with high
technology, practicality
and safety.

*comercial name

Water steam, ozone
steam* and aromatherapy

Skin hygiene and
aromatherapy with high
technology, practicality
and safety.

*comercial name

BEAUTY STEAM MAXX



Revolutionary design, which combines beauty and practicality in operation;

Accompanies cluster of LEDs for Chromotherapy
and Photobiomodulation, allowing the generation
of different colors;

More than 50 pre-programmed protocols;

It operates on two mains voltages - bivolt
(110V-/220V-);

Digital control of the water level of the reservoir,
against overflow;

Intelligent heating system that monitors the
temperature of the boiler providing more safety
and agility at the beginning of the operation;

Stem with 360° rotatable emitting tip and silicone
grip area for greater operator comfort;

Microcontroller with digital parameter controls;

Compartment for inserting flavorings.

ACNE;
ALOPECIA;
AROMATHERAPY;
CHROMOTHERAPY;
HYDRATION;
REJUVENATION;
TISSUE REPAIR.

ANVISA - 80212489004



You better
everyday