



ELEVATE YOUR CLINIC TO THE HIGHEST STANDARD
OF EXCELLENCE WITH HTM EQUIPMENT.
NEW TECHNOLOGIES, EFFICIENCY, DESIGN AND
PROFESSIONAL PERFORMANCE.

(+55) 19 3808-7741
www.htmeletronica.com.br
comex@htmeletronica.com.br
comercial@htmeletronica.com.br

Av. Rio Nilo, 209 - Jardim Figueira
CEP: 13904-380 - Amparo - São Paulo, Brasil



We reserve the right to correct any graphical errors.
Merely illustrative images. Internal distribution.





HTM

You better
every day

We are a reference industry in the development and manufacture of electronic equipment aimed at the aesthetics and physiotherapy segment in Brazil. With more than 20 years in the market, our differentials are based on high technology, exclusive design and great efficiency of our products. All equipment is designed and developed in its own factory, with 100% national technology!

Our modern industrial structure has 5 thousand m², built under the strictest standards of production and environmental responsibility. We are constantly expanding, always recognizing and valuing the ethical commitment and well-being of our employees.



Innovating in technology and creating high quality equipment is our mission. In this path of overcoming, we continue to support and value our dealers, customers and partners with the commitment to offer the market the best and most effective products, integrating scientific knowledge with practicality.

Today, HTM Eletrônica is a leader in innovation, results and safety, parameters that will undoubtedly follow with us in the future.



It has 4 body applicators, 2 flat and 2 anatomical, practical and ergonomic, which allow perfect coupling in different body regions;

10.4" touch screen display with greater touch accuracy, high durability and all interactivity of the Pro-Performance system;

Provides an advanced experience in managing your equipment through IoT technology, using the HTM Match application;

It has the HTM Track, which makes it possible to track the equipment if enabled for this function;

Exclusive remote software self-management system;

Microcontrolled with fully digital parameter controls;

It has the option of 1 specific applicator (optional) intended for the treatment of the pelvic muscles;

Possibility of using up to 4 simultaneous applicators, with individual intensity adjustments;

Free mode that provides adjustment of customized parameters;

It has the smart energy system;

Master-control security key;

18-month warranty;

It comes with carrying case;

Voltage: 220V;

7 Teslas Total Power - True and real.

PELVIC FLOOR DYSFUNCTION;
MUSCLE FLACCIDITY;
MUSCULAR HYPERTROPHY;
URINARY INCONTINENCE;
REDUCTION OF MEASURES;
MUSCLE TONING.

ANVISA - 80212489007

EMPOWER



High intensity focused
electromagnetic field

The technology in
favor of body and beauty.

BEAUTY DERMO MAXX



Vacuum therapy,
cupping therapy,
pressotherapy and
phototherapy

Advanced multifunctional
vacuum therapy , cupping
therapy, pressotherapy and
phototherapy technology.

Full range;

Exclusive design with new digital
pressure sensor;

Vacuum LED applicator containing: 3 green LEDs (300mW),
3 red LEDs (300mW) and 3 infrared LEDs (300mW), in
continuous and pulsed modes;

It presents the possibility of individual or combined use of
the LEDs, with or without the use of vacuum;

New key for Cupping Therapy, with a special and selective
suction system for applying suction cups;

New digital pressure sensor with variable negative pressure
adjustment from -1.5kPa to -80kPa**;

Pressotherapy with 2 output channels, with different modes
of operation, insufflation and insufflation direction, being
CHANNEL A for legs or arms and CHANNEL B for the
abdomen;

Offers the best body, facial and capillary vacuum therapy
applicators* ;

INFO menu and QR Code with access to use and care tutorial
for better use of all equipment resources;

More than 90 pre-programmed protocols and possibility to
record new ones;

Language selection (Portuguese, English and Spanish);

5" screen touch display and medical standard plug for easy
fitting and high durability;

Allows you to choose between treatments: Vacuum LED,
Vacuum Therapy , Phototherapy, Pressotherapy, Cupping
Therapy , Diamond Peeling and Buttock Pump *.

LOCALIZED ADIPOSITY AND CELLULITIS;
LYMPHATIC DRAINAGE;
FIBROSIS;
ANALGESIA;
RELAXING MASSAGE;
VENOUS STASIS;
LYMPHEDEMA AND EDEMA;
ACNE, STRETCH MARKS AND FACIAL AGING;
ALOPECIA AND MELASMA;
TISSUE REPAIR;
PHOTODYNAMIC THERAPY.

ANVISA - 80212480024
*Some accessories for performing these techniques are optional;
**Equivalent to -600mmHg.

TECARE

Tecartherapy and radiofrequency

Challenge your senses with the new generation of Tecartherapy.



New Tecartherapy and Radiofrequency system;

Revolutionary "hands free" technique": transforming technology into a new sensory experience;

Exclusive Tecartherapy applicators specially developed for the perfect fit of the hands: ideal for the "hands on" technique";

Treatments in the aesthetic and rehabilitation area;

touch display 5" screen;

INFO menu with usage and care tutorial for better use of all equipment resources;

Intelligent programming with targeted protocols according to the therapeutic modality and treatment region;

Innovative Pro Performance System, thought about the user experience;

Exclusive applicator for intimate technique with remote activation and intensity adjustment;

Multifrequency Technology: 0.6MHz, 1.2MHz and 2.4MHz;

Multi Layer Treatment : automatic frequency variation, reaching different tissue depths;

Different emission modes: continuous and pulsed;

Different options of capacitive and resistive tips for Radiofrequency: monopolar, bipolar, tetrapolar, hexapolar and decapolar;

Tecartherapy tips of different sizes: 30,50,70 mm.;

New plug for easy fitting and high durability;

Full range.

SAGGING OF THE SKIN;
FACIAL AGING AND WRINKLES;
ACNE SCARS;
STRETCH MARKS;
CELLULITIS;
LOCALIZED ADIPOSITY;
ACNE VULGAR;
INTIMATE AESTHETICS;
INCONTINENCE
URINARY;
ANALGESIA;
REHABILITATION.

ANVISA - 80212480015

TENS and FES;
Allows battery recharge;
2 channels;
33 pre-programmed protocols;
Possibility to save individual protocols;
Easy access superior outputs;
Full range;
High definition and contrast colorful OLED display.

ANALGESIA;
FUNCTIONAL REHABILITATION;
MUSCLE STRENGTHENING;
TISSUE REPAIR.

ANVISA - 80212480027

FISIO STIM



Currents

Portable excitomotor current device.

Currents

Full generator of excitomotor currents.

Currents

Generator of excitomotor currents.

STIMULUS PHYSIO MAXX



Full range;
Rotating menu;
Electroacupuncture protocol;
Currents: TENS, FES, HighForce, Russian, MediumForce, SoftForce, Interferencial 2000Hz and 4000Hz, Galvanic, Microgalvanic, MENS;
4 independent output channels;
51 pre-programmed protocols;
12 therapeutic currents (polarized and biphasic).

ANALGESIA;
CELLULITIS;
CONCENTRATED ADIPOSITY;
MUSCLE STRENGTHENING;
PERMEATION OF ACTIVES;
PHYSICAL REHABILITATION;
TISSUE REPAIR.

ANVISA - 80212480017



TENS, FES, HighForce, Russian, MediumForce, SoftForce (250mA);
4 channels;
47 pre-programmed protocols;
Possibility to save individual protocols;
Easy access frontal outputs;
High durability medical standard plug.

ANALGESIA;
DRAINAGE;
FUNCTIONAL REHABILITATION;
MUSCULAR FLACCIDITY;
TISSUE REPAIR.

ANVISA - 80212480017

STIMULUS PHYSIO 4 CANAIS



2 channels for currents;
150 pre-programmed protocols;

Body ultrasound: 1 and 3 MHz, 24W of power and ERA of 8cm²;

Facial ultrasound: 5MHz, 9W of power and ERA of 3cm² (optional);

Currents: TENS, FES, HighForce (1000Hz), Microgalvanic, MENS, Interferential 2kHz and 4kHz, Galvanic, Russian Current (2500Hz), MediumForce (4000Hz), SoftForce (8000Hz) and Lipolysis.

ANALGESIA AND REHABILITATION;
CELLULITIS AND LOCALIZED ADIPOSITY;
STRETCH MARKS AND TISSUE REPAIR;
STRENGTHENING MUSCULAR;
REJUVENATION;
WATER RETENTION;
ASSET PERMEATION.

ANVISA - 8021280001

1 and 3 MHz, 21W of power and ERA of 7cm²;

24 pre-programmed protocols;

Intelligent programming of application time calculation (3MHz);

Connection for combination therapy.

CELLULITIS;
LOCALIZED ADIPOSITY;
WATER RETENTION;
TISSUE REPAIR;
ANALGESIA;
PHONOPHORESIS;
REHABILITATION.

ANVISA - 80212480001

SONIC COMPACT MAXX



Ultrasound and therapeutic currents

Innovative intelligent system for combination therapy.



Microcontrolled digital ultrasound

High technology in portable ultrasound.

SONIC COMPACT 1-3 MHZ

STIMULUS-R

TENS, FES, HighForce (1000Hz), Russian (2500Hz), MediumForce (4000Hz), SoftForce (8000Hz), Drainage, Lipolysis, Galvanic, Microgalvanic and MENS;

10 channels for biphasic currents and 2 channels for polarized currents;

46 pre-programmed protocols;

Easy access front exits;

Durable medical standard plug;

Current for electroacupuncture.

CELLULITIS;
LOCALIZED ADIPOSITY;
MUSCLE FLAVORITY;
CUTANEOUS PERMEATION OF ASSETS;
WATER RETENTION;
TISSUE REPAIR;
PHYSICAL REHABILITATION;
ANALGESIA.

ANVISA - 802124800111



Complete multiplatform electrostimulation.



FLUENCE MAXX is a Phototherapy equipment with the highest LED and LASER technology. It presents a wide range of resources for multitherapy and exclusive accessories, with 27 options: Laser Pens with higher power, LED Cluster (circular and linear), LED Cluster + Laser, Cluster LED MAXX, Cluster LED + Laser MAXX, Face LED Phototherapy Mask, Face LED Phototherapy Mask + Neck and Bracelet for ILIB Technique.

FLUENCE MAXX

LED and Laser

High technology in
LED and LASER
multitherapies.



Exclusive design;

5" touch screen display;

28 accessory options;

INFO menu with tutorial for to use and care for better use off all the equipment's resources;

Permit the individual choice or union of different light sources, coherent and non-coherent, increasing the spectrum of light spectrum and the range of treatments;

More than 1,700* pre-programed protocols and possibility to record new;

Topocospic tip in the Laser Pens allowing to locate points for acupuncture;

Exclusive HTM Smart ID system for 904nm pens;

Medical standard plug of easy fitting and high durability;

Pen and cluster shutter button;

Exclusive bracelet for transcutaneous ILIB technique (optional);

LASER STOP and emergency button;

Laser/LED test device;

Full Range.

ACNE;
ACUPUNCTURE;
ALOPECIA;
ANALGESIA;
CELLULITIS;
FAT LOCATED;
HYDRATION;
LYMPHATIC DRAINAGE;
LYMPHEDEMA;
MELASMA;
PHOTODYNAMIC;
THERAPY;
REJUVENATION;
TISSUE AND INTRAVAGINAL REPAIR;
TRANSCUTANEOUS ILIB.

ANVISA - 80212480028

**Specific for each applicator, depending on the number of accessories purchased.*

Individual or associated therapies;

Option to include Laser pen, LED cluster and LED+Laser cluster;

More than 1,600 protocols*;

Exclusive HTM Smart ID system**.

ACNE;
ACUPUNCTURE - TRANSCUTANEOUS ILIB;
ALOPECIA;
ANALGESIA AND PHYSICAL REHABILITATION;
LOCALIZED ADIPOSITY AND CELLULITIS;
HYDRATION;
LYMPHEDEMA AND LYMPHATIC DRAINAGE;
MELASMA AND REJUVENATION;
TISSUE REPAIR;
PHOTODYNAMIC THERAPY.

ANVISA - 80212480028

*Depending on the number of accessories purchased.

**For 904nm pens.



FLUENCE

LED and Laser

Innovative LED and Laser technology.

Currents

Electrostimulator portable for speech therapy.



TENS and FES;

2 channels;

12 pre-programmed protocols;

Specific accessories for intraoral and extraoral application;

Easy access upper exits;

High resolution and contrast color OLED display.

ANALGESIA;
MUSCLE STRENGTHENING.

ANVISA - 80212480027

FONO STIM

